




StudioAtm_01.jpg



StudioAtm_02.jpg



StudioAtm_03.jpg



StudioAtm_04.jpg



StudioAtm_05.jpg



StudioAtm_06.jpg



StudioAtm_07.jpg



StudioAtm_08.jpg



StudioAtm_09.jpg



StudioAtm_10.jpg



StudioAtm_11.jpg



StudioAtm_12.jpg



StudioAtm_13.jpg



StudioAtm_14.jpg



StudioAtm_15.jpg



StudioAtm_16.jpg



StudioAtm_17.jpg



StudioAtm_18.jpg



StudioAtm_19.jpg



StudioAtm_20.jpg



StudioAtm_21.jpg



StudioAtm_22.jpg



StudioAtm_23.jpg




StudioAtm_24.jpg

Filename description:

Light-Probe

DH107LX



Resolution: 4096 x 4096 pixel

DH107LP


Resolution: 2048 x 2048 pixel

DH107LT

Resolution: 512 x 512 pixel

Light Probe Blured

DH107LB




Resolution: 512 x 512 pixel

Spherical Map

DH107SX

DH107SN



Resolution: 6248 x 3124 pixel

Resolution: 2048 x 1024 pixel

DH107_XXL

Resolution: 10.000 x 5.000 pixel

Vertical Cross


DH107VC



Resolution: 4686 x 6248 pixel

Horizontal Cross

DH107HC



Resolution: 6248 x 4686 pixel